Gerät zur Iontophorese  
Device for Iontophoresis

# IontoBravo

Medizinprodukt / Medical Device

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# Content

Please read the Operators Manual carefully before using the device and pay attention to all references !

## Page

### I. Device for Tap Water Iontophoresis IontoBravo

Delivery & Accessories	31
Explanation of Symbols	32
Introduction - What is <b>IontoBravo</b>	32
Indications	33
Contraindications	33
Technical Data / Medical Device Reporting	34
Safety Adjustments	35
Warnings	36
Advice on Application of Electrodes	37
Accessories	38
Maintenance / Warranty	39

### II. Operating- and Display-Elements of IontoBravo

Installation of the <b>IontoBravo</b> tap water system	40-41
Operating Elements / Connection	42
Battery Compartment with LOCK-button	43
The Display	44
The Info-line	45
Real Time Clock Setup	46
Language Selection	47

### III. The Therapy with IontoBravo

Program 01 - Direct Current	48
Program 02 - Pulsed Current	49
The Tap Water Treatment	50
Important Instructions for Tap Water Treatment	51
Drug Iontophoresis	52

### IV. Control of Therapy 53

### V. Batteries and Battery-Charger 54

Trouble Shooting Guide	55
Waste Disposal / Literature / Copyright	56
Electromagnetic Compatibility	61-63

## Delivery & Accessories

### Included in the package:

REF: 010E-188: 1 unit **IontoBravo D+P (Direct Current & Pulsed Current)**, 1 stand, 2 electrode plates aluminium (300 x 120 mm), 2 terry bags, 1 lead wire to connect the electrodes, 2 electrodes aluminium for axillary use (50 x 50 mm), 2 sponge-bags, 2 foam rolls, 4 rechargeable batteries type AA, 1 battery-charger, 1 carrying case (that constitute 2 water trays) and 1 operators manual

REF: 010E-189: 1 unit **IontoBravo DC (Direct Current)**, 1 stand, 2 electrode plates aluminium (300 x 120 mm), 2 terry bags, 1 lead wire to connect the electrodes, 2 electrodes aluminium for axillary use (50 x 50 mm), 2 sponge-bags, 2 foam rolls, 4 rechargeable batteries type AA, 1 battery-charger, 1 carrying case (that constitute 2 water trays) and 1 operators manual

REF: 010E-190: 1 unit **IontoBravo PC (Pulsed Current)**, 1 stand, 2 electrode plates aluminium (300 x 120 mm), 2 terry bags, 1 lead wire to connect the electrodes, 2 electrodes aluminium for axillary use (50 x 50 mm), 2 sponge-bags, 2 foam rolls, 4 rechargeable batteries type AA, 1 battery-charger, 1 carrying case (that constitute 2 water trays) and 1 operators manual


**IontoBravo** shall only be used with original accessories.

The device model is displayed at all times below the pulse.

### Spare parts:

REF

- 011F-117: 4 rechargeable batteries type AA
- 011E-118: 1 battery charger
- 011E-191: 1 carrying case Ionto
- 011E-192: 2 electrode plates aluminium 300 x 120 mm
- 011E-193: 2 terry bags
- 011E-194: 1 Ionto lead wire
- 011E-195: 2 electrodes aluminium for axillary use 50 x 50 mm
- 011E-196: 2 sponge bags Ionto for axillary use 70 x 100 mm
- 011E-197: 2 foam rolls
- 011E-198: 1 Ionto stand
- 011E-199: 1 operators manual **IontoBravo**

### Optional set for drug-iontophoresis:

REF:

- 011E-114: 2 rubber electrodes 45 x 50 mm
- 011E-115: 2 sponge bags 60 x 70 mm
- 011F-116: 2 velcro straps
- 011E-120: 1 lead wire to connect the electrodes

## Explanation of Symbols

Attention  
Warning

Type BF applied parts

Serial-Number  
of the product

Keep dry

Batch code

Manufacturer

Date of  
manufacture

CE-mark: confirms  
compliance with  
European standards

Not suited for persons  
with cardiac pacemaker

Read the instructions

## Introduction - What is IontoBravo

**IontoBravo** is a state-of-the-art tap water iontophoresis device designed to treat excessive sweating of the hands, feet and armpits in accordance with the recommendations of the German Dermatological Society (DDG). Patients, doctors and therapists are the intended operators of this unit. The treatment is painless with low side effects and the device can be used without supervision of a therapist and without help of a second person. Children under the age of 11 shall not use it without adult supervision. Numerous scientific studies have shown that this type of therapy has been a highly effective alternative to medication or surgery for many years. Treatment with **IontoBravo** is done in tap water using a simple hand and foot bath. For this, water trays are used, in which weak DC or pulsed current is passed. The current shall only be increased until the patient feels a slight tingling sensation. The weak current is sufficient to reduce the excretion of the sweat glands, without damaging them. In the armpit area, the current is applied by using wet sponge bags, if necessary. Optimal results are achieved with direct current, sensitive patients may use the monophasic pulsed current which is also effective. Models 188 and 190 (see page 31) offer the pulsed current function.

With an initial daily treatment lasting approximately 15 minutes, sweat secretion should normalize within 3-4 weeks. However, this process is reversible, so 1-2 maintenance sessions per week are recommended.

**IontoBravo** can also be used for transdermal drug delivery by drug-iontophoresis.

## Indications

**IontoBravo** can be used for several kind of hyperhidrosis treatment such as palmar, plantar or axillary hyperhidrosis.

The device can also be used to facilitate transdermal delivery of drugs.

## Contraindications

### **IontoBravo may under no circumstances be used in the following cases:**

- By patients with electronic implants like pumps or a pace maker because functional disturbances can occur
- By patients with myocardial damage (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
- In presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide
- When operating machinery with a potential of danger
- By patients with epilepsy
- During pregnancy
- If severe local inflammation or thrombosis is present
- By patients with febrile illnesses or infectious diseases
- When using metal-containing contraceptive spirals
- In case of significant coagulation disorders

### **Please pay unconditional attention to the following fundamental rules:**

- Do not use **IontoBravo** in places with skin defects or malignant diseases.
- Do not stimulate transcerebral (via the head).
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery).
- Do not stimulate in the area of the Sinus knot (heart).
- Do not stimulate on the Bulbous (eye).
- Do not stimulate in the Larynx-Pharynx-Area (throat / mouth).
- Do not position the electrodes directly on the spine.

### **Take special care with the use of the stimulation device:**

- if after an operation the healing process can be interrupted by a muscle contraction
- after acute bleeding or bone fractures with a tendency to bleed
- during menstruation
- in case of desensitized skin with nerve damage
- if you suffer from eczema
- if the area of treatment is located close to osteosynthesis material (metal in the body)
- if patients cannot use the device properly, due to mental or physical limitations

## Technical Data

### IontoBravo DC / IontoBravo PC:

- 1-channel device for iontophoresis with 2 programs.
- Intensity: 0 to max. 30 mA into 500 ohm, 1000 ohm or 1500 ohm load +/-2%
- Power supply: 4 rechargeable (NiMH) batteries 1,5 V Type AA
- Voltage max.: 60 V +/-2%
- Dimensions: L 141,6 mm x W 74 mm x H 40,4 mm
- Weight: 275 g (incl. batteries)

### IontoBravo PC:

- Waveform: monophasic positive, rectangular
- Frequency: 9900 Hz +/-2%
- Pulse width: 50-90  $\mu$ s +/-10% (selectable)

Firmware version 3.1 or later

**Classification:** Internally powered, constant current, continuous operation, Type BF applied part, medical device class IIa, IP 22

**IontoBravo** is protected against solid foreign bodies with a diameter of  $\geq 12.5$ mm and against dripping water at an incline of up to 15° and is rated IP22.

With regular average use, the expected service life of **IontoBravo** is 5 years. Service life of accessories is one year. The storage time of the batteries is approximately 2 years and the charging life is about 500 cycles.

Max. environmental conditions for transportation: -25 to +70 degree centigrade, 0-93% atmospheric humidity. Recommended environmental conditions for storage: -25 to +40 degree centigrade, 0-93% atmospheric humidity. The device shall not to be stored in direct sunlight. Allow up to one hour for the device to return to normal operating temperature if stored in conditions below +5°C or above +40°C.

Environmental conditions for use: +5 to +40 degree centigrade, 0-93% atmospheric humidity, air pressure 700-1060 hPa. During use the surface temperature of the device can be 41°C or higher, but never exceeds 48°C. Please place the device in the supplied stand if the ambient room temperature exceeds +34°C to avoid unnecessary discomfort.

This device is quality controlled and fulfills all legal requirements for medical products. Therefore this device carries the CE-sign „CE 0123“. The serial number is on the rear of the unit.

**Manufacturer:** MTR+ Vertriebs GmbH,  
Kamenzer Damm 78, 12249 Berlin - [www.mtrplus.com](http://www.mtrplus.com)

### Medical Device Reporting:

If there is a serious deterioration in the state of your health when using the medical device, report this to your distributor or to the manufacturer and to the competent authorities of your country.

## Safety Adjustments

1. The +button become inactive after one minute running time of a program to avoid an unintentional increase of intensity. This lock is switched off by pressing the +button three-times repeatedly. Increasing the intensity is possible again for another minute.
2. The maximum possible intensity setting is limited in the programs, depending on the planned therapy area, in order to prevent damage due to excessive current intensities. In programs for tap water iontophoresis, the voltage (in V) is adjusted and the current intensity changes according to the bodies resistance. Therefore it may happen that the set voltage values are not reached because the intensity (in mA) has already reached its limit.

These intensity limits are set as follows:

Mode	DC max.	PC max.
Feet	25 mA	30 mA
Hands	15 mA	20 mA
Armpit	8 mA	10 mA
Drugs	5 mA	7 mA

3. High current at lower voltages indicates that the skin resistance of the patient is very low due to excessive sweating. The device will display OVERLOAD if the low resistance is causing high current. The device will continue at a limited voltage and may not be able to reach the set or maximum voltage setting, depending on how low the resistance is.
4. The circuit is constantly monitored during operation to detect an open circuit immediately.
5. The voltage is kept constant during operation in order to reduce the fence effect to a comfortable level which cannot be felt.
6. If the circuit is interrupted during operation (i.e. removing a hand or foot from the water bath), the voltage will drop immediately to 0. CIRCUIT will flash on the display to prompt the user to close the circuit. When this is done, the voltage will slowly increase to the previously set value at 1 V per second.
7. If the skin is too dry, the resistance increases in the water bath. This is registered by the device and the term PROTECT flashes in the display, indicating that the treatment should be resumed at a later date when sweating returns.

## Warnings

- **IontoBravo** is type BF equipment and may not be immersed in liquids. The device is protected against debris larger than 12.5 mm in diameter from entering the cabinet. It is splash resistant and rated IP22. Do not use it in rain or in environments with excessive moisture. For indoor use only. 
- **IontoBravo** may not be used in presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- If using rechargeable batteries, be sure to use a CE approved battery charger. Never connect the **IontoBravo** directly to the mains, to a battery charger or to any other mains powered equipment.
- Burn marks could occur underneath the electrodes if used with surgical high frequency apparatus.
- The application of the electrodes for axillary treatment can increase the risk of cardiac fibrillation.
- If device operation occurs close to short wave, microwave, cellular or other wireless equipment, instability of the values of the unit could occur. The required minimum distance is 2 meters.
- Patients with metal and/or electronic implants may only use the unit with permission of an experienced doctor.
- Treatment in the areas of tumors is allowed only by instruction of the specialized physician.
- The **IontoBravo** may not be operated if the battery compartment is not closed. Before opening the battery compartment, make sure that there are no lead wires in the jacks at the head of the unit to eliminate any possibility of current flow.
- The operator should not touch the battery contacts & patient simultaneously.
- If you have swallowed any part of the device or any of its accessories, please contact a doctor immediately.
- **IontoBravo** should be used by children under the guidance of a doctor or therapist only.
- Drug iontophoresis must be done according to the instructions of a physician or clinician only.
- The use of gel-based electrodes could cause skin damage during drug iontophoresis.
- Metallic jewelry (rings, watches, bracelets etc.) have to be removed before treatment, since it can cause light burn marks.
- Minor skin lesions should be covered with vaseline or a greasy creme to isolate them from current flow.
- If the skin is particularly dry after treatment, if redness or even small cracks are observed, please apply a moisturizing lotion. If these complaints remain, please reduce the current intensity during the treatments. In heavier cases we recommend that you stop the therapy completely and consult your doctor.
- The treatment electrodes must always be covered by the supplied terry bags. Direct skin contact with the metal plates must be avoided.
- Axillary electrodes are allowed for use in wet sponge bags only.
- Two devices must not be used simultaneously on one patient.
- Don't use drug iontophoresis treatment in conjunction with ultrasound.

**Please, keep this device and its accessories stored out of the reach of children. Do not leave the lead wires in reach of children, since strangulation is possible due to the length of the lead wires.**

## Advice on Application of Electrodes

- Make sure, that the **IontoBravo** is switched off (no display visible on the LCD-Display).
- Connect the single poled ends of the lead wires to the electrodes plugs.
- The other end of the electrode lead wire will be plugged into the output jack A of the **IontoBravo**. Make sure, that all plugs fit completely to the electrode sockets as well as the device socket.
- Make sure that the protective film has been removed from the metal electrodes.
- The large metal electrodes for use with the water trays must always be completely covered by the terry bags.
- Moisten the sponge bags for axillary treatment carefully and do not compress. Then put the small metal electrodes into the bigger sponge bags (7 x 10 cm), so that the entire metal area is covered.
- The application area of the sponge bags intended for drug iontophoresis is  $6 \times 7 = 42 \text{ cm}^2$ , an intensity of 4 mA is recommended. If pulsed current is used with electrodes of this size, an intensity of up to 7 mA can be selected, depending on the selected pulse width. Please never use smaller electrodes with the **IontoBravo**.
- After each treatment please take the electrodes out of the sponge- or terry-bags and disconnect them from the lead wire. Then clean the electrodes as well as the bags with clear water.
- Finally allow the electrodes, terry bags and sponge bags to air dry.

### In general:

- Electrodes may not be altered.
- Do not pull the electrodes on the electrode lead wire.
- Apply the electrodes to clean and healthy skin only.
- If you are not using water baths, the hair should be shaved or shortened with scissors in case of strong bodily hair growth.
- For tap water iontophoresis please do not use ointments or creams in the area of treatment since these substances may have an isolating effect. Areas with small skin defects, that shall consciously be isolated, are excluded from this.
- Effective current densities of more than  $2 \text{ mA/cm}^2$  require an increased attention of the user with regard to the intensity adjustment.
- All materials that adhere to the skin from the sponge bags, electrodes or water can be removed with water and soap.
- Should serious skin irritation occur whilst using the **IontoBravo**, please interrupt the application and inform your doctor.

## Accessories

### Sponge- and Terry-bags

The sponge bags are made from commercial viscose and consist of a mixture of cellulose and cotton. The material is kind to the skin and meets the requirements of DIN EN 13432 (compostable material).

Before first use of any bag, please rinse with warm water once.

The material provides a high water absorption. Please make sure that the treatment is done with a wet bag only. If the sponges are slightly smaller and stiffer after drying, this does not effect the quality. Once they are soaked for therapy again, they get back their original flexibility.

If necessary, the sponge- or terry-bags can be cleaned with clean, warm water up to 60°. Please avoid the use of detergents or disinfectants, ironing, bleaching, chemical cleaning and tumble drying.

Residues of detergents or disinfectants could have

detrimental effects with monophasic pulses, due to ion transport.

**Please airdry all bags after use completely, don't put them back into the plastic bag when damp!**

The sponge- and terry-bags are intended for single patient use only.

Defect bags may not be used anymore and need to be replaced.

### Treatment tubs

The two plastic tubs of the carrying case can be separated from each other and serve as trays for hand or foot treatment. These trays can be cleaned with normal household cleaning products and are intended for use by a single patient only. Before being used again as a transport case, both trays should be thoroughly dried.

### Foam rolls

The two foam rolls can be placed on the short side of the treatment tubs to protect the wrists and forearms from uncomfortable pressure points. The material absorbs only little water and can be used for at least one year if used as intended. These rolls should be dried after use and are intended for use by a single patient only.

### Metal electrodes

The plate electrodes and the axillary electrodes are made of aluminium and contain no nickel components. If allergies still occur, please contact your doctor. Please use only clean water without any additives for cleaning. All electrodes are for single patient use only.

**Please remove the protective film before first use.**

## Maintenance

**IontoBravo** has been designed to be maintenance free. Nevertheless, please adhere to the following advice:

If needed, please clean **IontoBravo** - as well as the treatment tubs - with a mild detergent, water and a damp cloth. Please do not use alcohol based cleaning products. Do not expose **IontoBravo** to extreme humidity and do not submerge under water. A clean cloth, moistened with an alcohol free disinfection liquid can be used to clean the device. Please follow all instructions of the manufacturer of the liquid.

You may clean the lead wires with a damp cloth. You may also from time to time treat the lead wires with talcum powder in order to prevent cracks and to improve lead wire life.

The **IontoBravo** does not require calibration.

The **IontoBravo** is not repairable and contains no user serviceable parts.

No modification to the **IontoBravo** is allowed.

Remove the batteries from the **IontoBravo** if it is not used for extended periods of time.

After termination of the therapy, **IontoBravo** can be re-used by other patients. Please clean it before with a mild disinfectant cleaner.

## Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date). In case of failure, the **IontoBravo** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.

Excluded from this warranty are failures of the **IontoBravo** which result from excessive use, careless or improper handling as well as damage by force.



The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

The manufacturer reserves the right of technical changes, which may serve to increase treatment success, the functionality of the device or its accessories.

## Installation of the Tap Water System

Open the carrying case and take out the complete content. The two sides of the case can be separated - if necessary - and serve as treatment tubs for foot or hand therapy.

### Install the system as follows.

Connect the lead wire to the connectors of the electrodes. Please take care that the electrode connectors are covered completely.

Put the metal plates into the terry bags completely and take care that the connectors are near the opening of the bag.

Then lay down the terry bags with electrodes on the bottom of the water trays.



Connect the lead wire to the jack in the head of the device.

Put the unit into the stand on the opposite side of the body.

For hand treatment fix the foam rolls as wrist protection to the small side of the trays near your body.


Fill the trays with max. 1 cm clean and warm tap water without any additives.

We recommend having a bucket ready to pour the water out of the tubs after treatment.

## Installation of the Tap Water System



## Operating Elements

The TIME button allows to select the treatment time and can also be used to pause the program. In this case, it starts again after the break with a continuous increase up to the previously set intensity.

After reading the statistics, they can be deleted by holding the TIME-button.

## Connection

Jack to connect the electrode-lead wire

## Battery Compartment with LOCK-button

**IontoBravo** allows the doctor or therapist to lock one program for the patient to use. In this way the medical practitioner can be quite sure that the patient only uses the device with the correct program. The LOCK-button can be pressed with a pin of the lead wire or any other suited thin object. Pressing will fix the selected program so that different programs cannot be used anymore.


UNLOCK-button  
removes the  
program lock and  
shows global  
statistics

Scheme for placement of 4 batteries type AA

## The Display

The screen of **IontoBravo** consists of two different parts. The first part - above and below - includes fixed signs and symbols. All important informations for operation and therapy-monitoring are indicated here by numbers and letters. The second - shaded - area in the middle provides an optical impression of the waveform that is currently used in the particular program and names the used kind of **IontoBravo**.

## The Info-Line

The info-line in the upper part of the screen shows the following information:

### **Start information:**

IONTO DP = **IontoBravo** with Direct Current and Pulsed Current

IONTO DC = **IontoBravo** with Direct Current

IONTO PC = **IontoBravo** with Pulsed Current

### **Setup real time clock:**

SET HR - SET MIN - SET DAY - SET MTH - SET YR

### **Treatment information:**

DC FEET = Treatment of the feet with direct current

DC HANDS = Treatment of the hands with direct current

DC ARMS = Treatment of the armpits with direct current

DC DRUGS = Transdermal delivery of drugs by direct current

PC FEET = Treatment of the feet with pulsed current

PC HANDS = Treatment of the hands with pulsed current

PC ARMS = Treatment of the armpits with pulsed current

PC DRUGS = Transdermal delivery of drugs by pulsed current

CIRCUIT = The circuit of current is not closed. Maybe there is a defect of the lead wire - or one foot/hand is not in the tray.

PRG END = Treatment finished (end of program)

PAUSE = Interruption of program

OVERLOAD = Max. current limit is reached, due to very low skin resistance

PROTECT = Overtreatment protection - in case of very dry skin

More details in the info line are presented in the statistical analysis of treatments in chapter "Control of Therapy" (page 53).

## **Real Time Clock Setup**

Before first use of **IontoBravo** the “Real Time Clock” that is integrated should be set. It’s necessary for control of future therapy sessions. **IontoBravo** has a small battery that will support the clock with energy for several years. For this reason, there is no need to reset the clock, even if the 4 main batteries have been removed.

For setting time and date, please select any program 1-2 and press the **TIME**-button for at least 10 seconds. You’ll reach the configuration-mode. **HOUR** is indicated and you can set the hour with the **MODE**-button (00-23).

Press the left **+**button to confirm this choice and step forward to setting **MINUTE**. This is indicated on the screen as well, giving the opportunity to select the minutes via **MODE**-button in the same way as before.

This way you are able to set one after the other:

- **HOUR**
- **MINUTE**
- **DAY**
- **MONTH** and
- **YEAR**

Press the **PRG**-button on the frontside of the device to save the settings. All values will be automatically stored and the clock starts working.

## **Language Selection**

**IontoBravo** can be set to one of 6 different languages. If in the setting-mode of the RTC (see page 46) TIME is pressed for 5 seconds once again - now for the second time - the language selection screen is displayed.

The following options are available:

- German
- English
- French
- Spanish
- Italian
- Polish

Finally you store your choice with the **PRG**-button.

## Program 01

### **Treatment with direct current (P 01)**

This form of current may be slightly more effective for treatment, but may feel a little more uncomfortable for sensitive patients.

**Step 1:** The default treatment is set to 15 minutes. After selecting the program, the time display flashes and the treatment time can be individually determined via the TIME button between 5 and 30 minutes in steps of 5.

**Step 2:** After choosing the treatment time, the MODE button must be used to choose between foot-, hand- or armpit-treatment or drug iontophoresis. The selected mode is now shown in the infoline of the display.

**Step 3:** Now the treatment is started by pressing the left + button and the intensity is selected. In tap water iontophoresis modes this is done by selecting the voltage (V) displayed in the lower left corner of the display. The current intensity is calculated as a function of the skin resistance and is shown in the lower right corner of the display.

In the drug mode, select the intensity, which is displayed in the lower right corner (mA), directly and the voltage is calculated. The voltages determined in this way are displayed at the bottom left side.

Please always choose only enough volts or milliamps to feel the current. A higher current dosage is neither necessary nor useful. Both the voltage and the intensity have been factory limited in their maximum values to avoid any damage or unpleasant feeling due to excessive current intensities. As soon as one of the two limit values is reached, a further increase is not possible, so that the set values may not be reached.

The remaining time runs backwards towards 0 during the treatment on the display. If the treatment is interrupted at any time by pressing the PRG key, the device saves the remaining treatment time. The treatment can be continued by pressing TIME. Then the current slowly increases to the previously used intensity level. Pressing PRG again will end the therapy session.

## Program 02

### **Treatment with pulsed current (P 02)**

P 02 contains medium-frequency, monophasic rectangular pulses. This method is considered almost equivalent to DC. This option is recommended for more sensitive patients, especially children.

**Step 1:** After selecting the program, the pulse width (duration of each individual pulse) can be changed with the „ $\mu$ s“ key. Selectable are 50, 60, 70, 80 or 90  $\mu$ s. This corresponds approximately to a current flow of 50% to almost 90% of the treatment time. Due to the variable pulse width, the pulsed current can be designed according to individual ideas depending on the respective sensitivity. With a larger pulse width, a higher efficiency of the therapy is to be expected.

**Step 2:** The default treatment is set to 15 minutes. The treatment time can be individually determined via the TIME button between 5 and 30 minutes in steps of 5.

**Step 3:** After selecting the treatment time, the MODE button must be used to choose between foot, hand or armpit treatment or drug iontophoresis. The selected mode is now displayed in the infoline of the display.

**Step 4:** Now the treatment is started by pressing the left + button and the intensity is selected. In tap water iontophoresis modes this is done by selecting the voltage (V) displayed in the lower left corner of the display. The current intensity is calculated as a function of the skin resistance and is shown in the lower right corner of the display.

In the drug mode, select the intensity, which is displayed in the lower right corner (mA), directly and the voltage is calculated. The voltages determined in this way are displayed at the bottom left side.

Please always choose only enough volts or milliamps to feel the current. A higher current dosage is neither necessary nor useful. Both the voltage and the intensity have been factory limited in their maximum values to avoid any damage or unpleasant feeling due to excessive current intensities. As soon as one of the two limit values is reached, a further increase is not possible, so that the set values may not be reached.

The remaining time runs visibly backwards towards 0 during the treatment on the display. If the treatment is interrupted at any time by pressing the PRG key, the device saves the remaining treatment time. The treatment can be continued by pressing TIME. Then the current slowly increases to the previously used intensity level. Pressing PRG again will end the therapy session.

## The Tap Water Treatment

Install the **IontoBravo** Therapy System as described before (pages 40-41) on a dry, stable and level surface. For the hand or foot treatment fill the tubs with max. 1 cm of clean tap water, so that the terry bags are completely wet. Neither the back of the hand nor the back of the foot should be covered with water during the therapy. For armpit treatment, please use the included small metal electrodes in wet sponge bags.

**Be sure to remove all metal objects (watches, rings, bracelets, necklaces etc.) from your body.**

Now start **IontoBravo** with the central ON/OFF-button. The display will show the device type and its software version for one second. If it is a device with DC **and** pulsed current, the display then changes to the last used program. If you want to use a different current form, you can switch over with the PRG key. The treatment time is preset to 15 minutes, but can be changed with the TIME button. Then select the treatment mode (feet - hands - arms) with the right +/- buttons.

For **treatment of the hands** you select the current intensity with the left +/- buttons. Then immerse both hands concurrently with its entire area in the tubs. If the circuit has been closed in such a way, the current intensity will increase at a speed of max. 1 V per second until the set value is reached. It may take up to 2 minutes for the full intensity to be felt. Normally, however, significantly lower intensities are used. The current in the water bath is sufficient if you feel a slight tingling sensation. Excessive currents do not lead to better treatment success. For treatment of the hands, we recommend that the subsequent intensity adjustment should be carried out by a second person and only increase in increments of 1 V. If no second person is available, take your hand out of the water as completely as possible. In exceptional cases an uncomfortable, but harmless feeling can occur.

For **foot or axillary treatment**, pre-choice of intensity is possible, too. But we recommend, not setting the intensity until you have placed your feet in the tubs or the wet sponge bags in your armpits. The increase of intensity will happen at a maximum speed of 1 V per second. The intensity can be adjusted at any time during the treatment. But please note, that the safety features may require pressing the + button three times, if the +button has not been used for one minute or more.

During treatment you should not remove any part of the body from the tubs or avoid larger movements during the axillary treatment, as this may lead to unpleasant voltage peaks despite the protective circuit. When the treatment time has expired, the device automatically ramps down the current slowly and displays the end of the therapy session with PRG END in the display. Then you can take out your hands or feet from the tubs or remove the axillary electrodes.

## Important Instructions for Tap Water Treatment

Treatment observations have shown that the treatment success at the anode (+ Pol = red) is slightly better than at the cathode. If your problem is stronger on one side, please use the positive electrode indicated by the red connector on the lead wire for this side. In all other cases we recommend that you always change the polarity in the tubs or at the axillary electrodes between two treatment sessions. **Please never change during a therapy session.**

The current intensity while using the tubs is always sufficient for a good treatment result, when you feel a slight tingling sensation. You should not feel this tingling during therapy with axillary electrodes. If necessary, please adjust the intensity accordingly. With axillary electrodes a voltage below 5 V is normal.

If you need to take a part of your body out of the tub during the treatment, make sure you do it as complete as possible. If only small areas, such as toes or fingertips, are in contact with the water, voltage peaks can occur, which are harmless but may still be uncomfortable depending on the intensity.

If you want to interrupt the program, press the TIME button. The current intensity slowly drops to zero and the display will show PAUSE then. Press the TIME button again to resume the program and the power will slowly ramp up to the previously set level. **Never use the ON/OFF button during treatment**, since the sudden interruption of current may cause an uncomfortable, but harmless shock.

You may notice an increase in perspiration at the beginning of the therapy. This is normal and will improve after a few sessions.

Only use original accessories of MTR+, patient safety is not guaranteed when using third-party products.

Be careful not to allow the metal electrodes to come into direct contact with the skin during the therapy session. Therefore, always use the terry bags in the tubs and wet sponge bags during the axillary treatment.

Do not expose the device to excessive moisture or rain. If the device is exposed to excessive moisture or accidentally dropped in the water bath then turn it off, remove the battery cover and the batteries. Then dry everything thoroughly. Leave the device to air dry for 24 hours before use or storage. Ensure that all accessories are completely dry before storing them inside the carrying case.

Always pay attention to the battery charge status, which is displayed in the upper right corner of the display. If there is only one bar left, charge the batteries before you start treatment. If the state of charge drops down to one bar during a treatment, this therapy session can still be completed.


It is absolutely forbidden to mix batteries and rechargeable batteries. When replacing batteries, please always replace all 4 batteries.

## Drug Iontophoresis

**IontoBravo** is able to significantly improve the delivery of medication through the skin and allow deeper penetration into the tissues, so that a larger amount of the drug can be therapeutically effective. If it is noted on the leaflet of the medication that the drug is suitable for iontophoresis, it can be brought into the tissue with the aid of the device. In this case, the polarity of the substance must be named. A negatively charged drug is placed under the negative electrode (black connector), and vice versa a positive charged drug is placed under the positive electrode (red connector).

Only optional accessories from MTR+ may be used for this type of therapy (see page 31). The drug is applied to the skin and the appropriate rubber electrode is placed in a wet sponge bag on this area. The second electrode will be placed in a second wet sponge bag on the opposite side of the body. Both electrodes will be attached to the skin with a velcro strap or any tape and connected to the device via lead wire. Place the electrodes on healthy skin only and not over open wounds.

After switching on the **IontoBravo**, please select whether you want to work with direct current (P01) or pulsed current (P02). Then switch with the MODE button to DC DRUG. or PC DRUG. and start the therapy by pressing the + button, which also adjusts the intensity. The current should be about  $0.1 \text{ mA} / \text{cm}^2$  during a DC treatment. The application area of the sponge bags intended for drug iontophoresis is  $6 \times 7 = 42 \text{ cm}^2$ , so an intensity of 4 mA is recommended. If pulsed current is used with electrodes of this size, an intensity of up to 7 mA can be selected, depending on the selected pulse width. Please never use smaller electrodes with the **IontoBravo**.

Please thoroughly cleanse your skin before treatment and pay attention to unpleasant sensations during the treatment such as stinging or burning. If necessary reduce the intensity of the current or break off the treatment completely. Especially for drugs with a pain-reducing effect, carefully check your skin for redness or other irritation and visit your doctor as soon as possible if problems arise. Even if you do not feel the current, there are ions moving and it is no reason to increase the intensity. Scientific research has shown that a subliminal current dosage is advantageous in drug iontophoresis.

**Before you are using IontoBravo for drug iontophoresis you need to read the information that is delivered with the medication carefully.** Please check for the ingredients contained and follow all instructions regarding possible allergic reactions. If you are not sure, please ask your doctor before use.

## Control of Therapy

The **IontoBravo** records all treatment sessions and collects statistics. In this way doctors or therapists have the ability to better understand the possible reasons for the presence or the absence of treatment success. Longer sections of unintended use can be detected and prevented. The device stores a total of 60 days with 5 sessions per day. When 61 days of data are reached, consecutive days will overwrite the earliest days. The day number will continue increasing (61, 62, ...). For this reason the statistics will always show the actual data of the last 60 treatment days.

To read the memory, the UNLOCK-button in the battery-compartment (s.p. 43) must be pressed twice - if the device has been locked to a single program - or once, if the device has been unlocked. Now STATS. GL (for „global statistics“) will be displayed in the info line. If the device has not been used at all, the screen will show NO STATS.

### Global Statistics

Global statistics show the average values of all recorded sessions. The actual term is named in the info line and two small bars are flashing beside the respective value to make assignment easier. The change from one term to the next is done with the MODE-button. In this way the following information can be read one after the other:

- the most used PROGRAM
- the average used intensity in mA (AVG. MA)
- the number of treatment days (DAYS.USED)
- the number of treatment hours in total (HRS. USED)
- the average number of sessions per day (AVG. SESS)

By pressing the left +button (mA) the more detailed daily statistics are available. The screen shows STATS.DAY and the information regarding the first session of the first day can be seen. Pressing left +button (mA) = next day / Pressing left -button (mA) = previous day  
Pressing TIME = next session of the selected day  
Pressing right +/--buttons = scrolling through the session details.

### Session statistics

- number of treatment day (STATS.DAY)
- average intensity used in this session (AVG. MA)
- length of this session (SESS. LEN)
- used program (PROGRAM)
- date of this session (SESS.DATE)
- start time of this session (SESS.TIME)

Holding the left + or - button (mA) pressed for three seconds will return to „Global Statistics“. If the TIME-button is pressed for 5 seconds, DEL.STATS will be displayed and after confirmation with TIME - for 5 seconds once again - all statistics will be deleted and all individual programmed parameters are set to their default values.

## Batteries and Battery-Charger

**IontoBravo** is supplied with rechargeable batteries. If these batteries don't have sufficient charge to operate, you can recharge using the charger. For this purpose, please insert four rechargeable NiMH-batteries of type AA with right direction for polarity, as it is drawn in the charger and then plug it to the mains (230V). Charging starts immediately and after about 5 hours, the batteries are fully charged. There is no danger of overcharging because the device has a corresponding protection feature. While charging is in progress, the small yellow LED is blinking slowly. If this LED is continuously on, the batteries are charged completely. Now disconnect the unit from the power supply and then remove the charged batteries and put them into the battery compartment of **IontoBravo**.

If the yellow LED is flashing quickly, the charger stops charging because a wrong or damaged battery has been detected and needs to be removed.

If you want to use other types than the supplied batteries, please make sure that they have a minimum capacity of 1800 mAh or more - to avoid frequent charging.

### Safety Instructions:

- Don't charge batteries of different size (AA or AAA) at the same time.
- Please use NiMH-batteries only and do not try to charge normal batteries (primary batteries) - this could cause an explosion.
- Always charge 2 or 4 batteries, never 1 or 3 batteries. If you only load 2 batteries, please set them next to each other without a gap in the left or right charging slots, and not in the middle.
- Avoid the simultaneous charging of old and new batteries or charging of batteries with different capacity or of different manufacturers.
- Please never open batteries by force and don't throw them into a fire.
- The device may only be operated in enclosed, dry areas - it has to be protected from rain and moisture.
- It is normal that the batteries become hot when charging.
- If the power plug or the cabinet show any damage, the charger may not be operated.
- The charger may not be opened or altered, repairs shall be done by qualified persons only.
- Do not operate the charger unsupervised and disconnect it from the mains after use. Keep it out of reach of children or persons with reduced physical or mental capacity.
- Cleaning and maintenance may only be done after disconnection.

The address or website of the manufacturer can be found on the rear label of the charger. The actual charging current depends on the type and number of batteries to be charged. For loading the supplied four AA batteries it is 525 mA.

Please see also the instruction manual of the manufacturer.

## Trouble Shooting Guide

<b>Fault</b>	<b>Possible Cause</b>	<b>Solution</b>
CIRCUIT displayed on LCD during normal use.	Hand or foot not submerged in water.	Place hand or foot in water bath to complete the circuit.
	Lead wire not connected properly.	Check lead wire connections and ensure they are fully inserted.
	Lead wire is faulty.	Replace lead wire.
	Sponge bag is too dry.	Moisten sponge bag and repeat treatment.
	Device is damaged.	Return to the distributor for repair.
Device does not power on.	Flat battery or no batteries installed in device.	Check / Replace batteries.
	Battery installed in wrong orientation.	Remove batteries and install according to diagrams in battery compartment.
	Battery connector damaged or flattened.	Remove batteries, check metal parts in battery compartment for damage.
	Device is damaged.	Return to the distributor for repair..
Device powers off during use.	Batteries almost flat.	Check battery icon during use and replace batteries if low.
	Device is damaged.	Return to the distributor for repair.
OVERLOAD is displayed on LCD.	The device has detected that the patients skin has very low resistance due to very excessive sweating.	Current will be limited to prevent discomfort. Patient may continue treatment.
PROTECT is displayed on the LCD.	The device has correctly detected that the patients skin has high resistance due to very dry skin and does not require treatment.	Discontinue use until sweating returns. Test device on another person to confirm PROTECT is functioning correctly. Return device to manufacturer if PROTECT is displayed continuously even though patient is still sweating.
Treatment is painful (dull ache) or uncomfortable.	Intensity is too high for required treatment.	Reduce intensity (voltage) setting until treatment is comfortable. Note: Tingling sensations should be felt
Treatment is painful (sharp pain).	Patients skin has cuts or abrasions.	Discontinue use until skin is fully healed.
Skin becomes red after treatment.	Device is being used for periods which are too long for the patient.	Use device for a shorter treatment duration. Note: Some redness of skin is normal after treatment finishes.
Date is not recorded correctly in statistics.	RTC is not set correctly.	Repeat RTC-setup (see page 46)
	RTC battery for the internal clock is flat.	Return to manufacturer for servicing.

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In accordance with European regulations [1], used electrical and electronic equipment may not be discarded in household or unsorted municipal waste. Use separate collection facilities. Contact your local government or council for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the ground water and get into the food chain damaging the health and well being of people and animals.

[1] European Community directive 2012/19/EU on waste electrical and electronic equipment.

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A list of literature regarding the topics of this manual can be found on our website „[www.mtrplus.com](http://www.mtrplus.com)“.

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Reserved for changes.

Printing or copying of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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## Electromagnetic Compatibility

If the device is used according to the instructions given in this manual, it corresponds to the safety requirements for medical electrical systems IEC 60601-1. **IontoBravo** was checked and is in accordance with the limits set for medical products in the standard IEC 60601-1-2:2014. These limits serve to the safety of the user regarding electromagnetic disturbances.

The **IontoBravo** is a Class B device and is suitable for home healthcare environments as well as typical clinical environments.

Guidance and manufacturer's declaration - electromagnetic emissions		
The <b>IontoBravo</b> unit is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>IontoBravo</b> should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The <b>IontoBravo</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The <b>IontoBravo</b> is suitable for use in all establishments, including domestic and medical establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Please adhere to these warnings when using the device:

Don't use **IontoBravo** near active HF surgical equipment or RF shielded rooms used for magnetic resonance imaging, where the intensity of EM disturbances is high.

Use of this equipment adjacent to, or placement on top of other equipment should be avoided in case the output of the **IontoBravo** is effected. If this is necessary, Both the **IontoBravo** and the other equipment should be monitored to verify that they are both operating normally.

Use of accessories and lead wires other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used not closer than 30 cm (12 inches) to any part of the **IontoBravo**, including lead wires specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

## Electromagnetic Compatibility

<b>Guidance and manufacturer's declaration - electromagnetic Immunity</b>			
The <b>IontoBravo</b> is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>IontoBravo</b> should assure that it is used in such an environment			
<b>Immunity Test</b>	<b>IEC 60601 Test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact  ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact  ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Not Applicable	<b>IontoBravo</b> is powered by batteries
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV. 100 kHz repetition frequency	Not Applicable	<b>IontoBravo</b> is powered by batteries
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phases 0% $U_T$ for 0.5 cycle  And at single phase 0° 0% $U_T$ for 1 cycle 70% $U_T$ for 25/30 cycles 0% $U_T$ for 250 cycles 0% $U_T$ for 300 cycles	Not Applicable	<b>IontoBravo</b> is powered by batteries
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC61000-4-6	150 KHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % AM at 1 kHz	150 KHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % AM at 1 kHz	Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	
NOTE 1: $U_T$ is the a.c. mains voltage prior to application of the test level. NOTE 2: At 80 MHz the higher frequency range applies. NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

## Electromagnetic Compatibility

### Guidance and manufacturer's declaration - electromagnetic Immunity

The **IontoBravo** is intended for use in the electromagnetic environment specified below. The customer or the user of the **IontoBravo** should assure that it is used in such an environment

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
	450	380 – 390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
	870						
	930						
	1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
	5240						
	5785						

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